

DENTAL CATALOGUE



Arkay pak instruments (pvt.) Ltd.

*Manufacturers and Exporters of High Class Surgical, Dental, Manicure & Pedicure
Implements & Scissors of all sorts.*

Arkay pak instruments (pvt.) Ltd.

Manufacturers and Exporters of High Class Surgical & Dental Instruments

Philosophy

We have been making high standard surgical instruments since the year 1990 with the help and supervision of German technicians and are continuously exporting our instruments to our clients in Germany.

We have reorganized our production systems and certified as ISO - 9001 : 2008, ISO - EN 13485 : 2003, & cGMP. We have also got CE Marking Certification. We produce high quality Surgical, Dental, Veterinary, Manicure / Pedicure instruments and Scissors of all kind under the supervision of highly technical staff. We use high quality steel and our instruments are fault free as they are passed through under every necessary process of quality control.

We have the policy to provide:

- (i) high quality instruments*
- (ii) competitive rates*
- (iii) fastest possible delivery to our valued clients.*

INTERNATIONAL
CERTIFICATION
ISO EN 13485 : 2003



cGMP

ISO-9001-2008
APPROVED
CE MARK
Approved




CONTENTS

Description	Page #
Extracting Forceps	01-36
Root Elevators	37-59
Scalers, Probes & Spatulas	60-74
Oral & Maxillo-facial Surgery Instruments	75-95
Orthodontic Instruments, Orthodontic Pliers, Cutters, Impression Trays & Articulators	96-112
Index	113

Extracting Forceps

- English Pattern
- American Pattern



FINEST QUALITY DENTAL INSTRUMENTS

EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)

Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)




D-12-0001
1
 obere Schneide- und Eckzähne
 upper laterals and canines
 incisives et canines supérieures
 incisivos y caninos superiores




D-12-0002
2
 obere Schneide- und Eckzähne
 upper laterals and canines
 incisives et canines supérieures
 incisivos y caninos superiores




D-12-0003
3
 obere Schneide- und Eckzähne
 upper laterals and canines
 incisives et canines supérieures
 incisivos y caninos superiores



D-12-0004
4
 untere Schneide- und Eckzähne
 lower laterals and canines
 incisives et canines inférieures
 incisivos y caninos inferiores



D-12-0007
7
 obere Prämolaren
 upper premolars
 prémolaires supérieures
 premolares superiores



D-12-0008
8
 untere Prämolaren
 lower premolars
 prémolaires inférieures
 premolares inferiores




D-12-0013
13
 untere Prämolaren
 lower premolars
 prémolaires inférieures
 premolares inferiores

EXTRACTING FORCEPS (English Pattern)


Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)

Daviers (Forme Anglaise)
Forceps de Extraction (Forme Inglesa)




D-12-00135
135
 untere Prämolaren, für Kinder
 lower premolars, for children
 Premolaires inférieures, pour enfants
 premolares inferiores para niños

D-12-0017
17
 obere Molaren, rechts
 upper molars, right
 molaires supérieures droites
 molares superiores derechos



D-12-0018
18
 obere Molaren, links
 upper molars, left
 molaires supérieures gauches
 molares superiores izquierdos



D-12-0018A
18A
 obere Molaren
 upper molars
 molaires supérieures
 molares superiores



D-12-0018B
18B
 obere Molaren
 upper molars
 molaires supérieures
 molares superiores




D-12-0019
19
 obere dritte molaren
 upper third molars
 troisièmes molaires supérieures
 terceros molares superiores

EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)

Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)




EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)


Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)




ID-12-0029
29
 obere wurzeln
 upper roots
 racines supérieures
 raíces superiores




D-12-0029S
29S
 obere wurzeln, für kinder
 upper roots, for children
 racines supérieures, pour enfants
 raíces superiores, para niños




ID-12-0030
30
 obere wurzeln
 upper roots
 racines supérieures
 raíces superiores



D-12-0031
31
 untere wurzeln
 lower roots
 racines inférieures
 raíces inferiores



ID-12-0032
32
 untere wurzeln
 lower roots
 racines inférieures
 raíces inferiores



ID-12-0030S
30S
 obere wurzeln, für kinder
 upper roots, for children
 racines supérieures, pour enfants
 raíces superiores, para niños


EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)


Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)




D-12-0033
33
 untere wurzeln
 lower roots
 racines inférieures
 raíces inferiores



D-12-0033A
33A
 untere wurzeln
 lower roots
 racines inférieures
 raíces inferiores



D-12-0033B
33B
 untere wurzeln
 lower roots
 racines inférieures
 raíces inferiores



D-12-0033L
33L
 untere wurzeln
 lower roots
 racines inférieures
 raíces inferiores



D-12-0037
37
 oberer Schneidezähne, für Kinder
 upper incisors, for children
 incisives supérieures, pour enfants
 incisivos superiores, para niños




ID-12-0033S
33S
 untere Wurzeln, für Kinder
 lower roots, for children
 racines inférieures, pour enfants
 raíces inferiores, para niños


EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)


Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)




D-12-0038
38
 untere Schneide- und Eckzähne, für Kinder
 lower incisors and canines, for children
 incisives et canines inférieures, pour enfants
 incisivos y caninos inferiores, para niños




D-12-0039
39
 obere Molaren, für Kinder
 upper molars for children
 molaires supérieures, pour enfants
 molares superiores, para niños




D-12-0039R
39 R
 obere Molaren rechts, für Kinder
 upper molars right for children
 molaires supérieures, droites, pour enfants
 molares superiores, derechos, para niños



D-12-0039A
39 A
 obere Molaren, für Kinder
 upper molars for children
 molaires supérieures, pour enfants
 molares superiores, para niños



D-12-0040
40
 untere Molaren, für Kinder
 lower molars for children
 molaires inférieures, pour enfants
 molares inferiores, para niños



D-12-0039L
39 L
 obere Molaren links, für Kinder
 upper molars left for children
 molaires supérieures, gauches, pour enfants
 molares superiores, izquierdos, para niños

EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)


Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)



D-12-0041
41
 obere Wurzeln
 upper roots
 racines supérieures
 raíces superiores



D-12-0044
44
 obere Wurzeln
 upper roots
 racines supérieures
 raíces superiores



D-12-0048
48
 untere Wurzeln
 lower roots
 racines inférieures
 raíces inferiores



D-12-0045
45
 untere Wurzeln
 lower roots
 racines inférieures
 raíces inferiores



D-12-0049
49
 obere Wurzeln
 upper roots
 racines supérieures
 raíces superiores



D-12-0050
50
 obere Wurzeln
 upper roots
 racines supérieures
 raíces superiores

EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)

Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)



D-12-0051
51
 obere Wurzeln
 upper roots
 racines supérieures
 raíces superiores




D-12-0051A
51 A
 obere Wurzeln
 upper roots
 racines supérieures
 raíces superiores




D-12-0051L
51 L
 obere Wurzeln
 upper roots
 racines supérieures
 raíces superiores



D-12-0052
52
 obere Wurzeln
 upper roots
 racines supérieures
 raíces superiores



D-12-0052A
52 A
 obere Zähne und Resektion
 upper teeth and resection
 dents supérieures et résection
 dientes superiores y resección




D-12-0051S
51 S
 obere Wurzeln für Kinder
 upper roots, for children
 racines supérieures, pour enfants
 raíces superiores, para niños

EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)


Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)




D-12-0054
54
 Spaltzange für obere Molaren
 for separating upper molars
 à séparer molaires supérieures
 para separar molares superiores



D-12-0055
55
 Spaltzange für obere Molaren
 for separating upper molars
 à séparer molaires supérieures
 para separar molares superiores




D-12-0056
56
 Spaltzange für untere Molaren
 for separating lower molars
 à séparer molaires inférieures
 para separar molares inferiores



D-12-0065R
65 R
 obere Molaren, rechts
 upper molars, right
 molaires supérieures droites
 molares superiores derechos



D-12-0065L
65 L
 obere Molaren, links
 upper molars, left
 molaires supérieures gauches
 molares superiores izquierdos



D-12-0067
67
 obere dritte Molaren
 upper third molars
 troisièmes molaires supérieures
 terceros molares superiores



D-12-0067A
67 A
 obere dritte Molaren
 upper third molars
 troisièmes molaires supérieures
 terceros molares superiores


EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)

Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)




D-12-0067½R
67½ R
obere dritte Molaren, rechts
upper third molars, right
troisième molaires supérieures, droites
terceros molares superiores, derechos




D-12-0067½L
67½ L
obere dritte Molaren, links
upper third molars, left
troisième molaires supérieures, gauches
terceros molares superiores, izquierdos




D-12-0068
68
untere Wurzeln
lower roots
racines inférieures
raíces inferiores




D-12-0068A
68 A
untere Wurzeln
lower roots
racines inférieures
raíces inferiores



D-12-0073
73
untere Molaren
lower molars
molaires inférieures
molares inferiores



D-12-0073R
73 R
untere Molaren, rechts
lower molars, right
molaires inférieures, droites
molares inferiores, derechos



D-12-0073L
73 L
untere Molaren, links
lower molars, left
molaires inférieures, gauches
molares inferiores, izquierdos


EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)


Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)




D-12-0073A
73 A
 untere Molaren
 lower molars
 molaires inférieures
 molares inferiores




D-12-0073S
73 S
 untere Molaren
 lower molars
 molaires inférieures
 molares inferiores



D-12-0073B
73 B
 untere Molaren
 für Kinder
 lower molars
 for children
 molaires inférieures,
 pour enfants
 molares inferiores,
 para niños



D-12-0074
74
 untere Wurzeln
 lower roots
 racines inférieures
 raíces inferiores



D-12-0074D
74 D
 untere Wurzeln
 lower roots
 racines inférieures
 raíces inferiores



D-12-0073K
73 K
 untere Molaren
 für Kinder
 lower molars
 for children
 molaires inférieures,
 pour enfants
 molares inferiores,
 para niños


EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)


Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)




D-12-0074M
74M
 untere wurzeln
 lower roots
 racines inférieures
 raíces inferiores




D-12-0074N
74N
 untere wurzeln
 lower roots
 racines inférieures
 raíces inferiores



D-12-0074XN
74XN
 untere Wurzeln und eng-
 stehende schneidezähne
 lower roots and crowded
 incisors
 racines inférieures et
 incisives surdents
 raíces inferiores e incisivos
 supernumerarios



D-12-0075
75
 untere Prämolaren
 lower premolars
 prémolaires inférieures
 premolares inferiores



D-12-0076
76
 obere wurzeln
 upper roots
 racines supérieures
 raíces superiores




D-12-0076N
76N
 obere wurzeln
 upper roots
 racines supérieures
 raíces superiores

EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)

Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)




D-12-0076S
76S
 obere wurzeln
 upper roots
 racines supérieures
 raíces superiores



D-12-0079
79
 untere dritte Molaren
 lower third molars
 troisièmes molaires inférieures
 terceros molares inferiores



D-12-0079A
79A
 untere dritte Molaren
 lower third molars
 troisièmes molaires inférieures
 terceros molares inferiores



D-12-0081
81
 untere prämolaren und wurzeln
 lower premolars and roots
 Premolaires et racines inférieures
 premolares y raíces inferiores



D-12-0080
80
 Molarenwurzeln oben
 upper molar roots
 racines de molaires supérieure
 raíces de molares superiores

EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)

Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)



D-12-0086A
86A
 untere Molaren
 lower molars
 molaires inférieures
 molares inferiores




D-12-0086B
86B
 untere Molaren
 lower molars
 molaires inférieures
 molares inferiores




D-12-0086C
86C
 untere Molaren
 lower molars
 molaires inférieures
 molares inferiores




D-12-0089
89
 obere molaren, rechts
 upper molars, right
 molaires supérieures droites
 molares superiores derechos



D-12-0087
87
 untere Molaren
 lower molars
 molaires inférieures
 molares inferiores




D-12-0090
90
 obere molaren, links
 upper molars, left
 molaires supérieures gauches
 molares superiores izquierdos



EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)

Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)



EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)

Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)



D-12-0099 1/2
 # 99 1/2
 untere Molaren
 lower molars
 molaires inférieures
 molares inferiores



D-12-0099 1/2
 # 99 1/2
 untere Molaren
 lower molars
 molaires inférieures
 molares inferiores



D-12-0100
 # 100
 obere Wurzeln
 upper roots
 racines supérieures
 raíces superiores




D-12-0101
 # 101
 obere Pramolaren und wurzeln
 upper premolars and roots
 prémolaires et racines supérieures
 premolares y raíces superiores

EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)

Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)



D-12-0104
104
 obere Wurzeln
 upper roots
 racines supérieures
 raíces superiores




D-12-0102
102
 obere Molaren, rechts
 upper molars, right
 Molaires supérieures droites
 Molares superiores derechos




D-12-0103
103
 obere Molaren, links
 upper molars, left
 Molaires supérieures gauches
 Molares superiores izquierdos



D-12-0105
105
 untere Schneidezähne
 lower incisors
 incisiver inférieures
 incisivos inferiores



D-12-0106
106
 untere Schneidezähne
 lower incisors
 incisiver inférieures
 incisivos inferiores




D-12-0107
107
 obere Eckzähne
 upper canines
 canines supérieures
 caninos superiores

EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)


Daviers (Forme Anglaise)
Forceps de Extraction (Forme Inglesa)




D-12-0112
 # 112
 obere Wurzeln
 upper roots
 racines supérieures
 raíces superiores




D-12-0113
 # 113
 obere Wurzeln
 upper roots
 racines supérieures
 raíces superiores



D-12-0136
 # 136
 obere Prämolaren und Wurzeln
 upper premolars and roots
 prémolaires et racines supérieures
 premolares y raíces superiores



D-12-0137
 # 137
 untere Schneidezähne
 und wurzeln
 lower incisors and roots
 incisives et racines
 inférieures
 incisivos y raíces
 inferiores



D-12-0114
 # 114
 untere Wurzeln
 lower roots
 racines inférieures
 raíces inferiores




D-12-0143
 # 143
 untere Wurzeln
 lower roots
 racines inférieures
 raíces inferiores


EXTRACTING FORCEPS (English Pattern)

Extracting Forceps (English Pattern)
Zahnzangen (Englische Form)


Daviers (Forme Anglaise)
Forceps de Extracción (Forme Inglesa)




D-12-MD1
MD1
MEAD
 obere Frontzähne, Wurzeln
 upper anteriors, roots
 incisivos, caninos,
 racines supérieures
 incisivos, caninos,
 raíces superiores



D-12-MD2
MD2
MEAD
 obere Molaren
 upper Molars
 molaires supérieures
 molares superiores



D-12-MD3
MD3
MEAD
 untere wurzeln
 lower roots
 racines inférieures
 raíces inferiores



D-12-MD4
MD4
MEAD
 untere Molaren
 lower Molars
 molaires inférieures
 molares inferiores




D-12-0121
121
FELSCH
 obere dritte Molaren
 upper third molars
 troisièmes molaires supérieu
 terceros molares superiores



WISDOM EXTRACTING FORCEPS

Wisdom Extracting Forceps
Weisheitszahnzangen


Daviers pour dents de sagesse
Forceps para cordales



D-12-22 ½ L
Routurier
(# 22 ½ L)
Molaren und Weisheits-
zähne, unten, links
lower molars and
wisdoms, left
molaires et dents de sagesse
inférieures, gauches
molares y cordales inferiores,
lado izquierdo



ID-12-33 ½ L
Routurier
(# 33 ½ L)
Wurzeln unterer Molaren
und Weisheitszähne, links
lower molar and
wisdom roots, left
racines des molaires et dents
de sagesse inférieures, gauches
raíces de molares y cordales
inferiores, izquierdos



ID-12-33 ½ R
Routurier
(# 33 ½ R)
Wurzeln unterer Molaren
und Weisheitszähne, rechts
lower molar and
wisdom roots, right
racines des molaires et dents
de sagesse inférieures, droites
raíces de molares y cordales
inferiores, derechos

FORCEPS FOR CHILDREN

Forceps for children
Zahnzangen für Kinder

Daviers pour enfants
Forceps para niños



D-12-001S
Fig. 1 Baby
obere Schneidezähne
upper incisors
incisives supérieures
incisivos superiores



D-12-002S
Fig. 2 Baby
obere Prämolaren
upper premolars
prémolaires supérieures
premolares superiores



D-12-003S
Fig. 3 Baby
obere Molaren
upper molars
molaires supérieures
molares superiores




D-12-004S
Fig. 4 Baby
obere Wurzeln
upper roots
racines supérieures
raíces superiores



D-12-005S
Fig. 5 Baby
untere Schneidezähne
lower incisors
incisives inférieures
incisivos inferiores



D-12-006S
Fig. 6 Baby
untere Molaren
lower molars
molaires inférieures
molares inferiores




D-12-007S
Fig. 7 Baby
untere Wurzeln
lower roots
racines inférieures
raíces inferiores

ROOTS SPLINTER FORCEPS

Root Splinter Forceps
Wurzelsplitterzangen

Pinces Fragments
Pinzas Para Fragmentos



D-12-002W
Witzel
obere Wurzeln
upper roots
racines supérieures
raíces superiores



D-12-001W
Witzel
untere Wurzeln
lower roots
racines inférieures
raíces inferiores




D-12-001F
Ralk
Wurzelfragmente
root fragments
fragments de racines
fragmentos de raíces



D-12-002F
Stieglitz
Wurzelfragmente
root fragments
fragments de racines
fragmentos de raíces



D-12-003F
Stieglitz
Wurzelfragmente
root fragments
fragments de racines
fragmentos de raíces




D-12-004F
Arthur
Wurzelfragmente
root fragments
fragments de racines
fragmentos de raíces

EXTRACTING FORCEPS (American Pattern)

Extracting Forceps (American Pattern)
Zahnzangen (Amerikanische Form)

Daviers (Forme américaine)
Forceps de Extracción (Forme americana)



C -12A-0001
1
 obere schneidezähne und wurzeln
 upper incisors and roots
 incisives et racines supérieures
 incisivos y raíces superiores

D 12A-0002
2
 obere und untere wurzeln
 upper and lower roots
 racines supérieures et inférieures
 raíces superiores y inferiores

D -12A-00W1
W1
 winter
 obere schneide- und Eckzähne
 upper centrals and canines
 incisives et canines supérieures
 incisivos y caninos superiores

D -12A-0004
4
 obere und untere prämolaren
 upper and lower premolars
 prémolaires supérieures et inférieures
 premolares superiores y inferiores


D -12A-0005
5
 physick
 untere dritte Molaren
 lower third molars
 troisièmes molaires inférieures
 terceros molares inferiores

D -12A-0006
6
 spaltzange
 separating forceps
 pince à séparer
 forceps para separ

EXTRACTING FORCEPS (American Pattern)

Extracting Forceps (American Pattern)
Zahnzangen (Amerikanische Form)

Daviers (Forme américaine)
Forceps de Extracción (Forme americana)



D-12A-0008
8
 obere und untere dritte Molaren
 upper and lower third molars
 troisièmes molaires supérieures
 et inférieures
 terceros molares superiores y inferiores

D-12A-0013
13
 obere Schneidezähne
 upper incisors
 incisives supérieures
 incisivos superiores

D-12A-0007
7
 obere Wurzeln
 upper roots
 racines supérieures
 raíces superiores

D-12A-0010S
10S
 obere Wurzeln
 upper roots
 racines supérieures
 raíces superiores

D-12A-0010H
10 H
 obere Molaren
 upper molars
 molaires supérieures
 molares superiores

D-12A-0014
14
 untere Schneidezähne
 lower incisors
 incisives inférieures
 incisivos inferiores

EXTRACTING FORCEPS (American Pattern)

Extracting Forceps (American Pattern)
Zahnzangen (Amerikanische Form)

Daviers (Forme américaine)
Forceps de Extracción (Forme americana)



D-12A-0016
16
untere Molaren
lower molars
molaires inférieures
molares inferiores



D-12A-0015
15
Harris
untere Molaren
lower molars
molaires inférieures
molares inferiores



D-12A-0016S
16S
untere molaren, für kinder
lower molars for children
molaires inférieures, pour enfants
molares inferiores, para niños



D-12A-0018R
18R
Harris
obere Molaren, rechts
upper molars, right
molaires supérieures droites
molares superiores derechos



D-12A-0017
17
untere Molaren
lower molars
molaires inférieures
molares inferiores




D-12A-0018L
18 L
Harris
obere Molaren, links
upper molars, left
molaires supérieures gau
molares superiores izquie

EXTRACTING FORCEPS (American Pattern)

Extracting Forceps (American Pattern)
Zahnzangen (Amerikanische Form)


Daviers (Forme américaine)
Forceps de Extracción (Forme americana)




EXTRACTING FORCEPS (American Pattern)

Extracting Forceps (American Pattern)
Zahnzangen (Amerikanische Form)


Daviers (Forme américaine)
Forceps de Extraction (Forme americana)




D-12A-0030S
30S
obere zähne, für kinder
upper teeth, for children
dents supérieures, pour enfants
dientes superiores, para niños




D-12A-0032
32 *Parmly*
obere Prämolaren und Molaren
upper premolars and molars
prémolaires et molaires supérieures
premolares y molares superiores




D-12A-0032A
32A *Parmly*
obere Prämolaren und Wurzeln
upper premolars and roots
prémolaires et racines supérieures
premolares y raíces superiores



D-12A-0034
34
untere zähne, und wurzeln
lower anteriors, and roots
dents et racines, inférieures
dientes y raíces inferiores



D-12A-0035
35
obere Wurzeln
upper roots
racines supérieures
raíces superiores




D-12A-0036
36
untere zähne, und wurzel
lower anteriors, and roots
dents et racines, inférieures
dientes y raíces inferiores

EXTRACTING FORCEPS (American Pattern)

Extracting Forceps (American Pattern)
Zahnzangen (Amerikanische Form)


Daviers (Forme américaine)
Forceps de Extracción (Forme americana)




D-12A-0037
37
 untere Zähne und Wurzeln
 lower anteriors and roots
 dents et racines inférieures
 dientes y raíces inferiores




D-12A-0053R
53 R
 obere Molaren, rechts
 upper molars, right
 molaires supérieures droites
 molares superiores derechos



D-12A-0053 L
53 L
 obere Molaren, links
 upper molars, left
 molaires supérieures gauches
 molares superiores izquierdos



D-12A-0038
38
 obere Zähne und Wurzeln
 upper anteriors and roots
 dents et racines supérieures
 dientes y raíces superiores



D-12A-0062
62
 Schneidezähne, Prämolaren, Milchzähne
 incisors, premolars, deciduous teeth
 incisives, prémolaires, dents de lait
 incisivos, premolares, dientes de leche




D-12A-0065
65
 obere Schneidezähne und Wurzeln
 upper incisors and roots
 incisives et racines supérieures
 incisivos y raíces superiores


EXTRACTING FORCEPS (American Pattern)

Extracting Forceps (American Pattern)
Zahnzangen (Amerikanische Form)


Daviers (Forme américaine)
Forceps de Extracción (Forme americana)




D-12A-0069
69 Tomas
 obere und untere Wurzeln
 upper and lower roots
 racines supérieures et inférieures
 raíces superiores y inferiores




D-12A-0085
85 Tomas
 untere Eckzähne, Pramolaren, Molaren
 lower canines, premolars and molars
 canines, prémolaires, molaires infer.
 caninos, premolares, molares infer.



D-12A-0088R
88R Nevius
 obere Molaren, rechts
 upper molars, right
 molaires supérieures droites
 molares superiores derechos



D-12A-0085A
85A Tomas
 untere Eckzähne, Pramolaren, Molaren
 lower canines, premolars and molars
 canines, prémolaires, molaires infer.
 caninos, premolares, molares infer.




D-12A-0088L
88L Nevius
 obere Molaren, links
 upper molars, left
 molaires supérieures gauches
 molares superiores izquierdos


EXTRACTING FORCEPS (American Pattern)

Extracting Forceps (American Pattern)
Zahnzangen (Amerikanische Form)


Daviers (Forme américaine)
Forceps de Extracción (Forme americana)



D-12A-0099A
99A
Kells
 obere Schneide- u. Eckzähne, Prämolaren
 upper incisors, canines, premolars
 incisives, canines, prémolaires supér.
 incisivos, caninos, premolares super.




D-12A-0099C
99C
Kells
 obere Schneide- u. Eckzähne, Prämolaren
 upper incisors, canines, premolars
 incisives, canines, prémolaires supér.
 Incisivos, caninos, premolares super.



D-12A-0104
104
 obere dritte Molaren
 upper third molars
 troisièmes molaires supérieures
 terceros molares superiores



D-12A-0101
101
Hull
 Milchzähne
 Deciduous teeth
 dents temporaires
 dientes temporales




D-12A-00103
103
 untere Schneidezähne, Prämolaren
 lower incisors, premolars
 incisives prémolaires, racines infér.
 Incisivos premolares, raíces infer.


EXTRACTING FORCEPS (American Pattern)

Extracting Forceps (American Pattern)
Zahnzangen (Amerikanische Form)


Daviers (Forme américaine)
Forceps de Extracción (Forme americana)




D-12A-0150
150
Cryer
 obere Schneidezähne, Prämolaren
 und Wurzeln
 upper incisors, premolars, roots
 incisives, prémolaires, racines supér.
 Incisivos, premolares, raíces super.




D-12A-0150A
150A
Cryer
 obere Schneidezähne, Prämolaren
 und Wurzeln
 upper incisors, premolars, roots
 incisives, prémolaires, racines supér.
 Incisivos, premolares, raíces super.




D-12A-0150S
150S
Cryer
 obere Zähne, für Kinder
 upper teeth, for children
 dents supérieures, pour enfants
 dientes superiores, para niños



D-12A-0151
151
Cryer
 untere Zähne, für Kinder
 lower teeth, for children
 dents supérieures, pour enfants
 dientes superiores, para niños



D-12A-0151A
151A
Cryer
 untere Schneidezähne, Prämol
 und Wurzeln
 lower incisors, premolars, roots
 incisives, prémolaires, racines i
 Incisivos, premolares, raíces in




D-12A-0151S
151S
Cryer
 untere Schneidezähne, Prämol
 und Wurzeln
 lower incisors, premolars, roots
 incisives, prémolaires, racines i
 Incisivos, premolares, raíces in

EXTRACTING FORCEPS (American Pattern)

Extracting Forceps (American Pattern)
Zahnzangen (Amerikanische Form)


Daviers (Forme américaine)
Forceps de Extracción (Forme americana)




EXTRACTING FORCEPS (American Pattern)

Extracting Forceps (American Pattern)
Zahnzangen (Amerikanische Form)

Daviers (Forme américaine)
Forceps de Extracción (Forme americana)



D-12A-0222
222
 untere dritte Molaren
 lower third molars
 troisièmes molaires inférieures
 terceros molares inferiores



D-12A-0286
286
 obere Schneidezähne und Wurzeln
 upper incisors and roots
 incisives et racines supérieures
 incisivos y raíces superiores



D-12A-0287
287 Swing
 untere Molaren
 lower molars
 molaires inférieures
 molares inferiores



D-12A-0288
288
 untere Eckzähne und Premolari
 Lower canines and premolars
 canines, prémolaires inférieures
 caninos, premolares inferiores

EXTRACTING FORCEPS (American Pattern)

Extracting Forceps (American Pattern)
Zahnzangen (Amerikanische Form)


Daviers (Forme américaine)
Forceps de Extraction (Forme americana)



D-12A-000C1
Cohen
universalzahnzange
universal forceps
davier universal
forceps universal

D-12A-0010SN
10 S Nevius
obere Eckzähne
upper canines
canines supérieures
caninos superiores

D-12A-001AH
1 A Henahan
obere Schneidezähne
upper incisors
incisives supérieures
incisivos superiores



D-12A-0010HN
10 H Nevius
obere Eckzähne
upper canines
canines supérieures
caninos superiores

D-12A-0004S
4 Standard
obere dritte Molaren
upper third molars
troisiemes molaires supérieur
terceros molares superiores


EXTRACTING FORCEPS (American Pattern)

Extracting Forceps (American Pattern)
Zahnzangen (Amerikanische Form)


Daviers (Forme américaine)
Forceps de Extracción (Forme americana)



D-12A-00FH
FH Woodward
 untere Molaren
 Lower molars
 molaires inférieures
 molares inferiores



D-12A-0089
89
 obere Molaren, rechts
 Upper molars, right
 molaires supérieures droites
 molares superiores derechos



D-12A-00FS
FS Woodward
 untere Molaren
 Lower molars
 molaires inférieures
 molares inferiores



D-12A-0090
90
 obere Molaren, links
 Upper molars, left
 molaires supérieures gauchi
 molares superiores izquierdas


EXTRACTING FORCEPS (American Pattern)

Extracting Forceps (American Pattern)
Zahnzangen (Amerikanische Form)

Daviers (Forme américaine)
Forceps de Exacción (Forme americana)



D-12A-0023S
23 S
 untere Molaren
 Lower molars
 molaires inférieures
 molares inferiores



D-12A-0017SK
17 SK
 untere Molaren
 Lower molars
 molaires inférieures
 molares inferiores



D-12A-0150SK
150 SK
 untere zähne und wurzeln
 Upper anteriors and roots
 dents et racines inférieures
 dientes y raíces inferiores



D-12A-0151SK
151 SK
 obere zähne und Wurzeln
 lower anteriors and roots
 dents et racines supérieures
 dientes y raíces superiores



Arkay pak instruments (pvt.) Ltd.

Manufacturers and Exporters of High Class Surgical, Dental, Manicure & Pedicure Implements & Scissors of all sorts.

Office:

Bara Pathar, Near Church, Christian Town, Muslim Street Sialkot Pakistan.

Tel: 92-52-4272513-15(3 Lines)

Fax: 92-52-4272513

E-mail: into@arkaysurgical.com Web: www.arkaysurgical.com

ISO 9001:2008



cGMP

CE MARK
Approved

