

DENTAL CATALOGUE



Arkay pak instruments (pvt.) Ltd.

*Manufacturers and Exporters of High Class Surgical, Dental, Manicure & Pedicure
Implements & Scissors of all sorts.*

Arkay pak instruments (pvt.) Ltd.

Manufacturers and Exporters of High Class Surgical & Dental Instruments

Philosophy

We have been making high standard surgical instruments since the year 1990 with the help and supervision of German technicians and are continuously exporting our instruments to our clients in Germany.

We have reorganized our production systems and certified as ISO - 9001 : 2008, ISO - EN 13485 : 2003, & cGMP. We have also got CE Marking Certification. We produce high quality Surgical, Dental, Veterinary, Manicure / Pedicure instruments and Scissors of all kind under the supervision of highly technical staff. We use high quality steel and our instruments are fault free as they are passed through under every necessary process of quality control.

We have the policy to provide:

- (i) high quality instruments*
- (ii) competitive rates*
- (iii) fastest possible delivery to our valued clients.*

INTERNATIONAL
CERTIFICATION
ISO EN 13485 : 2003



cGMP



ISO-9001-2008
APPROVED
CE MARK
Approved



CONTENTS

Description	Page #
Extracting Forceps	01-36
Root Elevators	37-59
Scalers, Probes & Spatulas	60-74
Oral & Maxillo-facial Surgery Instruments	75-95
Orthodontic Instruments, Orthodontic Pliers, Cutters, Impression Trays & Articulators	96-112
Index	113

Scalers, Probes & Spatulas



FINEST QUALITY DENTAL INSTRUMENTS

CALIBRATED PROBES

Calibrated Probes
Tasxhenmeßsonden

Sondas Calibradas
Sondes Graduées



Probes Single End




CALIBRATED PROBES

Calibrated Probes
Tasxhenmeßsonden


Sondas Calibradas
Sondes Graduées

Probes-double End

Long lasting, easy-to-read, color coded, metal probes



Explorers




SCALERS


Scalers
Zahnreiniger

Instruments à détartrer
Limpiadores de sarro


McCalls



Columbia



Jacquette




EXCAVATORS


Excavators
Excavatoren

Excavadores
Excavateurs


Root Canal Excavators



Root Canal



Root Canal Spreaders-stainless Steel




SPREADERS


*Spreaders
Spreizer*

*Extesires
Ecarteurs*

Root Canal Spreaders



Endodontic Spreaders




PROBER & EXPLORERS

Prober & Explorers
Zahnsonden

Sondas y Exploradores
Sondas et Explorateurs

Root Canal Pluggers



Root Canal Explorer



GRACEY-CURETTES

Gracey-Curettes
Gracey-Küretten

Gracey-curetas
Gracey-curtttes



D-16-0066
Gracey 1/4-short shank for
incisors and cuspids



D-16-0067
Gracey 1/4-short contra-angle shank for
incisors and cuspids



D-16-0068
Gracey 5/8-medium shank for
cuspids and molars



D-16-0069
Gracey 7/8-medium contra-angle shank for
bicuspid and molars, buccal and lingual surfaces



D-16-0070
Gracey 9/10-long contra-angled for molars
with more acute blade angle



D-16-0071
Gracey 11/12°-triple angulation for easy access to mesio-buccal
and mesial-lingual surfaces of bicuspid and molars



D-16-0072
Gracey 11/12°-triple angulation for easy access to
disto-posterior surfaces of bicuspid and molars



D-16-0073
Gracey 15/16-Angled shank, designed to reach
mesial surfaces of posterior teeth




D-16-0074
Gracey 15/16-Highly angled shank for reaching
distal surfaces of posterior teeth

UNIVERSAL CURETTES


Universal Curettes
Universal-Küretten

Curetas unersales
Curettes unverselles


204 Series



Sickles



Barnhart




FILLING INSTRUMENTS

Filling Instruments
Füllungsinstrumente

Instruments pour plombage
Instrumentos para obturaciones

Burnishers-double End















Burnishers-single End



BLACK'S CUTTING INSTRUMENTS

Black's Cutting Instruments
Schneidende instrmente nach black

Instrumentos cortantes de black
Instruments tranchants de black


	Shooshan A	D-16-0094
	Bredall 4	D-16-0095
	21	D-16-0096
	21B	D-16-0097
	Frahm H/H	D-16-0098
	Hollenback 3	D-16-0099
	Hollenback 3S	D-16-0100
	Walls 3	D-16-0101
	Ward 1	D-16-0102
	Ward 2	D-16-0103
	Caulk 5/7	D-16-0104
	Levy 7	D-16-0105

PLASTIC FILLING INSTRUMENTS & EXCAVATORS


Plastic Filling Instruments & Excavators
Instrumente für plastische füllungen, Excavatoren

Instrumentos para obturaciones plásticas, Excavadores
Instruments pour plombages plastiques, Excavateurs

Plastic Filling




Excavators



PLASTIC FILLING INSTRUMENTS

Plastic Filling Instruments
Instrumente für plastische füllungen












Instrumentos para obturaciones plásticas
Instruments pour plombages plastiques



WAX AND MODELLING CARVERS

Wax And Modelling Carvers
Wachs-und modellierinstrumente

Instruments á cire et á modeler
Talladores para cera y para modelar

	Black 89/92	D-16-0130
	Cleoid/Discoïd	D-16-0131
	3/6	D-16-0132
	5C	D-16-0133
	6C	D-16-0134
	7C	D-16-0135
	3T	D-16-0136
	IPC Interproximal	D-16-0137
	IPC-A (Angled)	D-16-0138
	4T	D-16-0139
	5T	D-16-0140









FILLING INSTRUMENTS

Filling Instruments
Füllungsinstrumente

Instruments pour plombage
Instrumentos para obturaciones

	Sherman Packer Smooth, Serrated	D-16-0141
	Balshi Packer Smooth, Serrated	D-16-0142
	Balshi Packer Angled	D-16-0143
	Round Angled Smooth, Serrated	D-16-0144
	Round Straight Smooth, Serrated	D-16-0145
	Taper Angled Smooth, Serrated	D-16-0146
	Taper Straight Smooth, Serrated	D-16-0147


Pluggers / Condensers

	Black 1/2 Smooth, Serrated	D-16-0148
	Hollenback 1 Smooth, Serrated	D-16-0149
	Mort 2 Smooth	D-16-0150
	Smith Smooth, Serrated	D-16-0151
	Detroit 1 Smooth	D-16-0152
	Ward 1 Smooth	D-16-0153
	Ward 2 Smooth	D-16-0154
	Ward 3 Smooth	D-16-0155

MIXING SPATULAS

*Mixing Spatula
Anmischspatel*

*Spatule à Malaxer
Espatula para mezclar*



Wax 7A
D-16-0156

Wax 7
D-16-0157

CMI 24
D-16-0158

36
D-16-0159

24
D-16-0160

24A
D-16-0161

22
D-16-0162



Arkay pak instruments (pvt.) Ltd.

Manufacturers and Exporters of High Class Surgical, Dental, Manicure & Pedicure Implements & Scissors of all sorts.

Office:

Bara Pathar, Near Church, Christian Town, Muslim Street Sialkot Pakistan.

Tel: 92-52-4272513-15(3 Lines)

Fax: 92-52-4272513

E-mail: into@arkaysurgical.com Web: www.arkaysurgical.com

ISO 9001:2008



cGMP

CE MARK
Approved

