

DENTAL CATALOGUE



Arkay pak instruments (pvt.) Ltd.

*Manufacturers and Exporters of High Class Surgical, Dental, Manicure & Pedicure
Implements & Scissors of all sorts.*

Arkay pak instruments (pvt.) Ltd.

Manufacturers and Exporters of High Class Surgical & Dental Instruments

Philosophy

We have been making high standard surgical instruments since the year 1990 with the help and supervision of German technicians and are continuously exporting our instruments to our clients in Germany.

We have reorganized our production systems and certified as ISO - 9001 : 2008, ISO - EN 13485 : 2003, & cGMP. We have also got CE Marking Certification. We produce high quality Surgical, Dental, Veterinary, Manicure / Pedicure instruments and Scissors of all kind under the supervision of highly technical staff. We use high quality steel and our instruments are fault free as they are passed through under every necessary process of quality control.

We have the policy to provide:

- (i) high quality instruments*
- (ii) competitive rates*
- (iii) fastest possible delivery to our valued clients.*

INTERNATIONAL
CERTIFICATION
ISO EN 13485 : 2003



cGMP

ISO-9001-2008
APPROVED
CE MARK
Approved



CONTENTS

Description	Page #
Extracting Forceps	01-36
Root Elevators	37-59
Scalers, Probes & Spatulas	60-74
Oral & Maxillo-facial Surgery Instruments	75-95
Orthodontic Instruments, Orthodontic Pliers, Cutters, Impression Trays & Articulators	96-112
Index	113

***Orthodontic Instruments,
Orthodontic Pliers, Cutters,
Impression Trays, Articulators***



FINEST QUALITY DENTAL INSTRUMENTS

ORTHODONTICS INSTRUMENTS

Orthodontics Instruments
Orthodontische Instrumente

Instrumentos Ortodonticos
Instruments Orthodontiques

BRACKET HOLDER TWEEZERS



D-20-0001
Bracket Placer
Tweezers
(Extra Fine Tip)



D-20-0002
Easy Access
Tweezers



D-20-0003
Slipper
Tweezers



D-20-0004
Direct bonding
bracket holders
135mm



D-20-0005
Direct bonding
bracket holders
136mm



D-20-0006
Direct bonding
bracket holders
135mm



D-20-0007
Deluxe Bracket
Holder



D-20-0008
Direct bonding
bracket holders
140mm



D-20-0009
Deluxe Bracket
Holder Tweezers



D-20-0010
Deluxe Bracket
Holder Tweezers



D-20-0011
Direct Bond
Bracket Holder
Tweezers



D-20-0012
Bracket
Tweezers
Slot Aligner



D-20-0013
Direct Bond
Bracket Holder
Tweezers, Small Tip



D-20-0014
Exact Bracket
Holder Tweezers


ORTHODONTICS INSTRUMENTS

Orthodontics Instruments
Orthodontische Instrumente

Instrumentos Ortodonticos
Instruments Orthodontiques



D-20-0015 Ortho Light Band Pusher Scaler



D-20-0016 Ortho Light double ended ligature director



D-20-0017 Ortho Light double ended ligature director scaler



D-20-0018 Ortho Light double ended ligature scaler/director
tivate coil spring are remove modules with one instruments




D-20-0019 Ortho band pusher scaler




D-20-0020 Ortho Light Band Pusher Scaler




D-20-0021 Ortho ligature director single ended



D-20-0022 Ortho Light double ended ligature director scaler



D-20-0023 Ortho Light double ended ligature scaler/director
tivate coil spring are remove modules with one instruments



D-20-0024 Ortho Light double ended ligature director



D-20-0025 Ortholite Ligature Director Scaler



D-20-0026 Double ended cement spatula



D-20-0027 Ortho light ligature director double ended




D-20-0028 Ortho light scaler single ended

ORTHODONTICS INSTRUMENTS

Orthodontics Instruments
Orthodontische Instrumente

Instrumentos Ortodonticos
Instruments Orthodontiques



D-20-0029 Ortho light scaler double ended




D-20-0030 Ortho light scaler double ended



D-20-0031 Ortho light scaler double ended



D-20-0032 Twiel-on elastic module instrument



D-20-0033 Twiel-on elastic module instrument



D-20-0034 Schure Band Seater



D-20-0035 Ortho light scaler



D-20-0036 Ortho light scaler



D-20-0037 Ortho light scaler heavy



D-20-0038 Ortho instrument




D-20-0039 Ortho instrument

ORTHODONTICS INSTRUMENTS

Orthodontics Instruments
Orthodontische Instrumente

Instrumentos Ortodonticos
Instruments Orthodontiques




D-20-0040 Band seater auto cleavable




D-20-0041 Band seater auto cleavable



D-20-0042 Ortho instrument



D-20-0043 Heavy duty band pusher



D-20-0044 Positive Positioner



D-20-0045 Polar Bear Ice Instruments



D-20-0046 Distal Bender

PLIERS FOR ORTHODONTICS & PROSTHETICS

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



D-20-0047
T.C
angled direct
bond bracket
remover



D-20-0048
T.C
Pin holding plier



D-20-0049
T.C
tweed loop
forming plier



D-20-0050
T.C
straight direct
bond bracket
remover



D-20-0051
T.C
flat nose plier



D-20-0052
bird beak plier

PLIERS FOR ORTHODONTICS & PROSTHETICS

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



D-20-0053
T.C
V bend stop plier



D-20-0054
T.C
Utility plier
weingart



D-20-0055
T.C
standard ligature
cutter



D-20-0056
T.C
V bend stop plier



D-20-0057
T.C
mini pin and
ligature cutter



D-20-0058
T.C
3-prong clasp
adjustable plier


PLIERS FOR ORTHODONTICS & PROSTHETICS

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



D-20-0059
T.C
Distal End Cutter
Safety hold is ideal for
cutting wires within
the patient's mouth.



D-20-0060
T.C
angled direct
bond bracket
remover



D-20-0061
T.C
Direct
bonding bracket
remover (wide) plier



D-20-0062
T.C
Hook crimping
pliers
angled straight



D-20-0063
T.C
Pin and ligature
cutter




D-20-0064
T.C
Direct bracket
remover



PLIERS FOR ORTHODONTICS & PROSTHETICS

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik



Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis





D-20-0065
 Aderer Clamp three jaw wire bending plier for contouring and bending wire up to .30"/0.76mm.




D-20-0066
 Bird beak plier for bending round, square and rectangle wire up to .30"/0.76mm.




D-20-0067
 Bird beak plier for light wire bending round, square and rectangular up to .016"/0.41mm.




D-20-0068
 Double Rounded jaw bending pliers, hooks and wire up to .020"/0.51mm.



D-20-0069
 Adams pliers for hooks and bending wire up to .030"/0.76mm.



D-20-0070
 Johnson pliers for crown band contouring



D-20-0071
 Falcon beak contouring plier.




D-20-0072
 For placement and removal of arch wires, pins and other auxiliaries.



D-20-0073
 For placement and removal of arch wires, pins and other auxiliaries in angle design



D-20-0074
 Bending Plier. In order to skirt molares bands. Non grooved.



D-20-0075
 contouring Plier.



D-20-0076
 Bending plier. In order to skirt molares bands. With grooved. Recommended wire .016", .018" and .022" 0.41mm, 0.46mm and 0.56mm


PLIERS FOR ORTHODONTICS & PROSTHETICS

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik



Pincès pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis




D-20-0077
 Nance loop closing plier, for bending of stop loops and adams' clasps. Wire up to .022" x .028", 0.7mm.



D-20-0078
 Loop forming in four length 3,4,5 & 6 mm wire capacity 0.22", 56mm.





D-20-0079
 Arch bending plier design for convenient handling of square or rectangular wire. .020" / 0.51mm.



D-20-0080
 In order to skirt bands and crowns.




D-20-0081
 Wires up to .020" / 0.7mm.




D-20-0082
 Hard wire cutter up to .022" x .028" .56mm x 71mm.




D-20-0083
 Round and concave plier for loop forming and banding wire up to .20" / .051mm.



D-20-0084
 Bending for round square or rectangular wire up to 0.30"/0.76mm



D-20-0085
 Plier for double bars.




D-20-0086
 for double back and triple back bands in 0.30" / .76mm & 0.36" / .91mm lingual arch wires.

PLIERS FOR ORTHODONTICS & PROSTHETICS

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pincas pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



D-20-0087
 Angle
 14.5 cm, 5 1/2"
 Bandformzange
 Band forming plier
 Pince à former les bagues
 Alicata para formar bandas




D-20-0088
 Peeso
 13 cm, 5"




D-20-0089
 Peeso
 13 cm, 5"

Kronenbänderzangen
 Collar Pliers
 Pincas pour bandes à couonnes
 Alicates para cuellos



D-20-0090
 Peeso
 14.5 cm, 5 1/2"

Kronenstreckzangen
 Crown stretching pliers
 Pincas à étendre les couonnes
 Alicates para ensanchar coronas




D-20-0091
 Peeso
 14.5 cm, 5 1/2"



D-20-0092
 Peeso
 14.5 cm, 5 1/2"

Zang flach/rund
 Plier flat and round
 Pince plate et ronde
 Alicata plano y redondo



D-20-0093
 16 cm, 6 1/4"

Universal-Kramponzangen
 Pin Roughening-Cutting and Bending Pliers
 Pincas universelle pour crampons
 Alicates para pone asperas



D-20-0094
 16 cm, 6 1/4"

PLIERS FOR ORTHODONTICS & PROSTHETICS

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



D-20-0095
13 cm, 5"



D-20-0096
Schwarz
13 cm, 5"




D-20-0097
Schwarz
13 cm, 5"




D-20-0098
Schwarz
14 cm, 5 1/4"


Draht- und Biegelbiegezangen
Arch and Spring Bending Pliers
Pinces à plier les fils et les barres
Alicates para doblar alambre y barras



D-20-0099
12.5 cm, 4 3/4"



D-20-0100
Angle
12.5 cm, 4 3/4"



D-20-0101
Adams
12.5 cm, 4 3/4"



D-20-0102
13 cm, 5"

Drahtbiegezangen
Wire Bending Pliers
Pinces à plier les fils
Alicates para doblar alambre

PLIERS FOR ORTHODONTICS & PROSTHETICS


Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



D-20-0103
 Young
 13 cm, 5"

Dahtbiegezange
 Wire Bending Plier
 Pince à courber les fils
 Alicate para doblar alambre



D-20-0104
 Tryfus
 16.5 cm, 6 1/2"

Kronenformzange
 Crown Forming Plier
 Pince pour former les couronnes
 Alicate para formar coronas




D-20-0105
 Johnson
 12 cm, 4 3/4"


Konturzangen
 Contouring Pliers
 Pince à buteroller
 Alicate para contornear



D-20-0106
 Reynolds
 13 cm, 5"




D-20-0107
 Gordon
 13.5 c, 5 1/4"



D-20-0108
 USA Modell
 13.5 cm, 5 1/4"

Konturzangen
 Contouring Pliers
 Pincas à buteroller
 Alicates para contornear



D-20-0109
 Goslee
 13.5 cm, 5 1/4"




D-20-0110
 Frevert
 13 cm, 5"

PLIERS FOR ORTHODONTICS & PROSTHETICS


Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis




D-20-0111
14 cm, 5 1/2"

Flachzange zum Einsetzen
von Einschraubspitzen
Flat nose plier for screw points
Pince plates pour pointes à vis
Alicate de picos chatos
para puntas a rosca




D-20-0112
16 cm, 6 1/4"


Flachzangen
Flat nose pliers
Pinces plates
Alicates de picos chatos



D-20-0113
16 cm, 6 1/4"




D-20-0114
14 cm, 5 1/2"




D-20-0115
14 cm, 5 1/2"

Flachzangen
Flat nose pliers
Pinces plates
Alicates de picos chatos



D-20-0116
13.5 cm, 5 1/2"



D-20-0117
13.5 cm, 5 1/2"

NEEDLE HOLDERS

*Nadelhalter
Porta-agujas*

*Porte-agujas
Portaghi*



D-20-0118
Mathieu Needle Holders
(Wide Tip) 14 cm



D-20-0119
Mathieu Ligature
Placing



D-20-0120
Mathieu Needle Holders
(Narrow Tip) 14 cm



D-20-0121
Mathieu Forceps
With Hook



D-20-0122
Mini Mathieu
Needle Holder



D-20-0123
Mathieu
Elastic Tip



D-20-0124
Mathieu Narrow Tip Double
Spring (T.C. Tip) Gold Handle






D-20-0125
Mathieu Medium Tip Double
Spring (T.C. Tip) Gold Handle

IMPRESSION TRAYS-STAINLESS STEEL

Impression Trays-Stainless Steel
Abdruckloffel S.S.


Cubetas para impresiones de acero inox
Porte-empreintes en acier inoxydable




IMPRESSION TRAYS-STAINLESS STEEL

Impression Trays-Stainless Steel
Abdruckkloffel S.S.


Cubetas para impresiones de acero inox
Porte-empreintes en acier inoxydable



D-20-0129
Impression Trays
Set of 8 Pieces
Perforated Stainless
steel Coe Type



D-20-0130
Impression Trays
Set of 8 Pieces
Plain Stainless Steel
Rim Lock Type for
Senior Persons



D-20-0131
Impression Trays
Set of 8 Pieces
Perforated Stainless
steel Coe Type
for Senior Persons

ARTICULATORS

Articulators Artikulations

Artocilar Articuler



D-20-0132
Stainless Steel
Finishing



D-20-0133
Golden Black Metallic



D-20-0134
Golden Orange Metallic



D-20-0135
Brown



D-20-0136
Royal Blue



D-20-0137
Silver Black Metallic



D-20-0138
Silver Blue Metallic



D-20-0139
Silver Green Metallic



D-20-0140
Copper Black Metallic



Arkay pak instruments (pvt.) Ltd.

Manufacturers and Exporters of High Class Surgical, Dental, Manicure & Pedicure Implements & Scissors of all sorts.

Office:

Bara Pathar, Near Church, Christian Town, Muslim Street Sialkot Pakistan.

Tel: 92-52-4272513-15(3 Lines)

Fax: 92-52-4272513

E-mail: into@arkaysurgical.com Web: www.arkaysurgical.com

ISO 9001:2008



cGMP

CE MARK
Approved

